



Scan me

ICONREHAB™

# whizziXbalance

FLIGHT  
SAFE  
BATTERIES

## Height adjustable, Auto Levelling & Folding Compact Powerchair

Are you looking for a high-performance, smart and foldable **electric wheelchair** combining comfort, control and safety? Look no further!

The **Whizzi X Balance Powerchair** is engineered for modern mobility – offering advanced electronic features, intelligent safety systems, app control function and exceptional usability in a lightweight and portable design

## Smart Technology & Intelligent Control

- **Electronic Postural Balance System:** Automatically adjusts and levels the seat on slopes or inclines to maintain upright posture and improve user stability.
- **One-Touch Electronic Folding & Unfolding:** Fold or unfold the powerchair with a single button press – perfect for travel and compact storage.
- **Adjustable Seat Height (Electronic):** Raise or lower the seat electronically for easier transfers, table access, or better visibility & inclusion.
- **20° Electrically Adjustable Backrest:** Recline at your convenience to reduce fatigue and enhance posture.



## Whizzi™ X Balance Features:

- Electrical Auto levelling Postural Balance Seat
- Electrically Adjustable Seat Height & 20° Recline
- Electrical Auto Fold Function
- Adjustable Leg Rest Support
- Adjustable Armrest Height
- 115Kg user weight capacity
- 17 km range between charges
- 2 x 220W motors
- 3 part Disassembly for easy transportation
- 82cm turning circle ideal for tight spaces
- iOS or Android Full Function Remote Control App



Breaks down into 3 main parts for easy transport



Freewheel  
Mode



Flight-safe  
battery option



2 x 220W power



15.5A Lithium  
battery



Puncture free  
tyres



One-touch  
electric fold



**HARAKA MOBILITY**  
JOY of MOVEMENT

## Haraka & Dubai Mobility

Unit 3, 416 Nadd Al Hamar, Dubai, United Arab Emirates.

Tel: [+971 4898 0300](tel:+97148980300) Whatsapp: [+971 55 116 6459](tel:+971551166459)

<https://dubaimobility.com>

**New!**

ICONREHAB™

# whizziXbalance

## Whizzi X Balance: Power, Portability and Peace of Mind

The Whizzi X Balance Electronic Wheelchair is more than just a mobility aid – it's your smart, stylish and reliable partner for everyday life. Whether at home, on the go, or exploring outdoors, it's built to give you the freedom you deserve.

### Advanced Technology Meets Everyday Freedom

Introducing the Whizzi X Balance Powerchair, a revolutionary folding electronic wheelchair engineered to redefine your mobility experience. Combining intelligent technology, ergonomic comfort and portable convenience, it empowers you to live independently, confidently and on your own terms.

### Effortless Stability and Exceptional Comfort on Every Journey

Conquer slopes and uneven surfaces with the Electronic Postural Balance Automatic Seat Levelling, designed to automatically adjust your seating angle when tackling inclines up to 15°. Stay upright, secure and balanced at every movement.

### Total Control, Enhanced Safety and Reliable Performance

Take charge with the fully featured hand controller and pressure-sensitive joystick, offering precise steering and intuitive speed control. With 5 speed levels, you can smoothly navigate indoor spaces, busy pavements and open areas.

### Personalise your Whizzi X Balance

Personalise your X Balance to reflect your personality with optional coloured rear shrouds

**Find out more or  
book your test drive**



Scan me



**HARAKA MOBILITY**  
JOY of MOVEMENT

### Haraka & Dubai Mobility

Unit 3, 416 Nadd Al Hamar, Dubai, United Arab Emirates.

Tel: [+971 4898 0300](tel:+97148980300) Whatsapp: [+971 55 116 6459](tel:+971551166459)

<https://dubaimobility.com>